

USER MANUAL

The MedSource Total Lifter is a transport unit rated to transport up to 450 kg. The unit must be placed so the handle support straps are to the side opposite the patient. The Total Lifter is for use by professional staff for lifting and transporting / transferring patients. It is recommended to use at least 4 persons to provide safe transport / transfer, thus reducing the possibility of accident or injury. Grab handles in line with the people lifting on opposite sides of the unit (DO NOT CROSS HANDLES). The Total Lifter is specifically for use when there is no need for spinal stabilization or traction of the patient being transported. If the Total Lifter is damaged, frayed, cut or soiled replace with new unit.

Directions

The MedSource Total Lifter is designed as an emergency manual transfer sheet to be used by a minimum of four people, with 6 people being recommended. The Total Lifter lifts, moves and extracts victims and patients.

Step 1: The sheet is positioned full length under the patient.

Step 2: Each person takes hold of two (one per hand) handles.

Step 3: Using only recommended / approved and proper safety technique the patient is lifted and moved accordingly.

Size & Weight

107 x 183 cm; 900 g

Material

Polyester Fibers - non-woven and woven with antimicrobial properties.

Weight Rating

Supports up to 450kg

Single handle load rating is 110 kg

Fluid Resistance

Resistant to the following fluids:

Oils, Alcohol, Mild Cleaning Fluids, Transmission Fluids and Brake Fluids

Packaging

1 per poly bag

10 per carton

	Consult instructions for Use
	Non sterile
	Caution, consult accompanying documents
	Manufacturer
	EEC Representative
	Catalogue number
	Do not reuse

Total Lifter

MS-95200

REF

 MedSource
8600 Shelby Court, Chanhassen, MN 55317 USA
952-472-0131
www.medsourcelabs.com