



MedSource
Respiratory & Airway

O2 Regulators



MS-010R25

MS-010R15

All Brass, 0-25 LPM

The MedSource 0-25LPM Oxygen Regulator is constructed of all brass internal and external parts. This unit provides flow rates of 0-25 LPM for "D" size oxygen cylinders and features two DISS ports. The "T" Style handle is fitted with a screw stop to prevent handle loss. Ignition tested according to CGA/ASTM standards.

Aluminum Body, 0-15 LPM

The MedSource 0-15 LPM Oxygen Regulator features an ergonomic "easy grip" dial and large easy-to-read numbers. The dial is constructed of impact resistant TPE plastic. The "T" style handle is fitted with a screw top to prevent handle loss. Ignition tested according to CGA/ASTM standards.

Rev. 1

www.gomedsource.com

Item	Description	UOM
MS-010R25	O2 Regulator, All Brass, 0-25 LPM	EA, 20 EA/CS
MS-010R15	Oxygen Regulator, Aluminum Body, 0-15 LPM	EA, 20 EA/CS