



MedSource
Diagnostics

Digital Thermometer



The MedSource digital thermometer is well equipped to give an accurate read through oral, underarm, or rectal usage. It has a 20mm x 7.5mm LCD Screen. Our thermometer also has a beeper function and automatic shut-off to conserve energy and prolong battery life. A 1.5v button battery and 5 disposable sheaths come included.

Rev. 1

www.gomedsource.com

Item #	Description	UOM
MS-121000	Digital Thermometer	EA, 12EA/BX 24 BX/CS 244 EA/CS