



**MedSource**  
Respiratory & Airway

# Reusable Silicone Resuscitator

The MedSource Reusable Silicone Resuscitators are complete with mask, reservoir bag and oxygen tubing.

## Features:

- Made from semi-transparent silicone
- Easy to disassemble & assemble for cleaning, disinfecting and sterilization at 134 degrees
- Equipped with pressure limitation valve for patient safety
- Textured surface for firm grip and effective ventilation

## Specifications:

- Maximum stroke: 1500ml
- Maximum resuscitator volume: 2000ml
- Oxygen reservoir volume: 2600ml
- Inspiratory resistance: 50 L/min
- Expiratory resistance: 5 L/min



Rev. 1

[www.gomedsource.com](http://www.gomedsource.com)

Item #	Description	UOM
MS-6110	Reusable Silicone Resuscitator (BVM), Adult	10 EA/CS
MS-6120	Reusable Silicone Resuscitator (BVM), Pediatric	10 EA/CS
MS-6130	Reusable Silicone Resuscitator (BVM), Infant	10 EA/CS