



MedSource

I.V. Infusion

I.V. Extension Sets

MS-83089		8", Luer Slip Lock, Slide Clamp. Approximate Priming Volume 2.0 mL	100 EA/CS, 6 CS/MC, 600 EA/MC
MS-83090		6", Needleless PRN, Split Septum Y-Site, Slide Clamp. Approximate Priming Volume 1.5 mL	100 EA/CS, 6 CS/MC, 600 EA/MC
MS-83091		6", Needleless Clear PRN, Roberts Clamp, Rotary Male Luer Lock. Approximate Priming Volume 2.0 mL	100 EA/CS, 6 CS/MC 600 EA/MC
MS-83092		8", Needleless PRN, Roberts Clamp. Approximate Priming Volume 2.0mL	100 EA/CS, 6 CS/MC, 600 EA/MC
MS-83094		10", Needleless PRN, Roberts Clamp Approximate Priming Volume 2.0 mL	100 EA/CS, 6 CS/MC, 600 EA/MC
MS-83097		19", Flow Regulator, Split Septum Y-Site, Slide Clamp, Male Luer Lock. Approximate Priming Volume 3.2 mL	48 EA/CS, 6 CS/MC, 288 EA/MC
MS-83216		21", Male Luer Lock. Approximate Priming Volume 2.3 mL	50 EA/CS, 6 CS/MC, 300 EA/MC
MS-83377		6", Needleless Clear PRN, Split Septum Y-Site, Rotary Male Luer Lock. Approximate Priming Volume 2.0 mL	50 EA/CS, 6 CS/MC, 300 EA/MC
MS-83378		8", Needleless Clear PRN, Roberts Clamp, Rotary Male Luer Lock. Approximate Priming Volume 2.0 mL	100 EA/CS, 6 CS/MC, 600 EA/MC

Our I.V. Extension Sets are available in multiple lengths to provide maximum flexibility when administering an I.V.

Rev. 1

www.gomedsource.com