



MedSource
Gear Bags

Deluxe Oxygen Bag



Deluxe O2 Bag Specification:

- Size: 27" x 12" x 10"
- Color: Red, Royal Blue

The Deluxe Oxygen Bag is padded and features two zippered pockets on the front of the bag with a reflective stripe for high visibility in low light emergency situations. There are also two zippered compartments on the sides for extra supply storage. The bag also features a removable tray under the lid with two hook-and-loop dividers providing for flexible storage configurations. The tray sits over the lower storage compartment facilitating easy access to a size "D" oxygen cylinder and supplies.

Rev. 1

www.gomedsource.com

Item #	Description	UOM
MS-B3312	Oxygen Bag, Deluxe, Red	EA, 4 EA/CS
MS-B3313	Oxygen Bag, Deluxe, Royal Blue	EA, 4 EA/CS