



**MedSource**  
I.V. Infusion

# Tourniquet



The MedSource tourniquet is a constricting or compressing device used to control venous and arterial circulation to an extremity for a period of time. Our Latex-Free Textured Tourniquets have a unique “grippy” surface to better hold the patient’s arm when drawing blood. Made from a blend of materials that closely resemble the physical properties of latex, without the worries associated with latex allergies. Our tourniquets measure 1” x 18” and are available rolled or unrolled.

Rev. 1

[www.gomedsource.com](http://www.gomedsource.com)

Item #	Description	UOM
MS-17301	Tourniquet, Latex-Free, 1' x 18", Disposable, flat	EA, 1000 EA/CS
MS-17302	Tourniquet, Latex-Free, 1' x 18", Disposable, rolled	EA, 250 EA/BG, 4 BG/CS, 1000 EA/CS