



**MedSource**  
Respiratory & Airway

# ET Tubes



The MedSource Endotracheal Tubes are ideal for oral intubation and airway management. Constructed of kink-resistant Medical Grade PVC, these ET Tubes also feature a radio-opaque black line. On cuffed versions, the pilot balloon allows for a tactile sense of cuff volume. Available cuffed or uncuffed, and with or without a pre-loaded stylet.

**Packaged Sterile, 10 EA/BX, 10 BX/CS, 100 EA/CS**

**This product is not made with natural rubber latex.**

**Cuffed:**

**Uncuffed:**

**Uncuffed, with Stylet:**

**Cuffed, with Stylet:**

Item #	Description	Item #	Description	Item #	Description	Item #	Description
MS-23250	5.0mm, Cuffed	MS-23120	2.0mm, Uncuffed	MS-23325	2.5mm, Uncuffed	MS-23450	5.0mm, Cuffed
MS-23255	5.5mm, Cuffed	MS-23125	2.5mm, Uncuffed	MS-23330	3.0mm, Uncuffed	MS-23455	5.5mm, Cuffed
MS-23260	6.0mm, Cuffed	MS-23130	3.0mm, Uncuffed	MS-23335	3.5mm, Uncuffed	MS-23460	6.0mm, Cuffed
MS-23265	6.5mm, Cuffed	MS-23135	3.5mm, Uncuffed	MS-23340	4.0mm, Uncuffed	MS-23465	6.5mm, Cuffed
MS-23270	7.0mm, Cuffed	MS-23140	4.0mm, Uncuffed	MS-23345	4.5mm, Uncuffed	MS-23470	7.0mm, Cuffed
MS-23275	7.5mm, Cuffed	MS-23145	4.5mm, Uncuffed	MS-23350	5.0mm, Uncuffed	MS-23475	7.5mm, Cuffed
MS-23280	8.0mm, Cuffed	MS-23150	5.0mm, Uncuffed	MS-23355	5.5mm, Uncuffed	MS-23480	8.0mm, Cuffed
MS-23285	8.5mm, Cuffed	MS-23155	5.5mm, Uncuffed			MS-23485	8.5mm, Cuffed
MS-23290	9.0mm, Cuffed	MS-23160	6.0mm, Uncuffed			MS-23490	9.0mm, Cuffed
MS-23295	9.5mm, Cuffed					MS-23410	10.0mm, Cuffed
MS-23210	10.0mm, Cuffed						