



MedSource
Respiratory & Airway

Nasal Cannula



The Nasal Cannula is designed for patients requiring medium or low concentrations of oxygen. The over the ear design maintains proper positioning of nasal tips, while allowing complete freedom of movement. The star-lumen supply tubing is durable and kink resistant. All varieties of nasal cannula are individually packaged in poly bags.

Oxygen Concentration Delivery: 24-44%
Recommended Flow Rate: 1-6 LPM

Rev. 1

www.gomedsource.com

Item #	Description	UOM
MS-24003	Nasal Cannula, Non-Flared Tip, 7'Tubing	EA, 50 EA/CS
MS-24004	Nasal Cannula, Flared Tip, 7'Tubing	EA, 50 EA/CS
MS-24005	Nasal Cannula, Non-Flared Tip, 14'Tubing	EA, 50 EA/CS
MS-24006	Nasal Cannula, Non-Flared Tip, 25'Tubing	EA, 25 EA/CS
MS-24007	Nasal Cannula, Non-Flared Tip, 50'Tubing	EA, 25 EA/CS
MS-24008	Nasal Cannula, Non-Flared Tip w/Univ. Con. , 7'Tubing	EA, 50 EA/CS
MS-24101	Nasal Cannula, Pediatric, 7'Tubing	EA, 50 EA/CS
MS-24203	Nasal Cannula, Adult, 7'Tubing w/Soft Tips	EA, 50 EA/CS
MS-24211	Nasal Cannula, Pediatric, 7'Tubing w/Soft Tips	EA, 50 EA/CS
MS-24212-U	Nasal Cannula, Infant , Soft Tips w/Univ. Con, 7'Tubing	EA, 50 EA/CS