



MedSource
Respiratory & Airway

Tracheostomy Mask



The MedSource Tracheostomy Mask allows for the delivery of oxygen therapy to patients who have had a tracheostomy.

The mask is constructed of molded transparent, Medical Grade PVC. It comes equipped with a swivel adapter and a spandex terylene adjustable strap for added comfort.

This product is not made with natural rubber latex.

Rev. 1

www.gomedsource.com

Item #	Description	UOM
MS-25100	Tracheostomy Mask, Adult	EA, 50 EA/CS