



MedSource
Respiratory & Airway

LED Standard Laryngoscope Blades

Our Disposable Laryngoscope blades feature LED bulbs that deliver a COLD, SAFE, and BRIGHT light source.

The Stainless Steel blade and aluminum heel block result in an instrument that is stronger than a plastic blade.



Rev. 1

www.gomedsource.com

Item #	Description	UOM
MS-464100	Disposable Standard Macintosh Blade, Size 00	EA
MS-46410	Disposable Standard Macintosh Blade, Size 0	EA
MS-46411	Disposable Standard Macintosh Blade, Size 1	EA
MS-46412	Disposable Standard Macintosh Blade, Size 2	EA
MS-46413	Disposable Standard Macintosh Blade, Size 3	EA
MS-46414	Disposable Standard Macintosh Blade, Size 4	EA
MS-464200	Disposable Standard Miller Blade, Size 00	EA
MS-46420	Disposable Standard Miller Blade, Size 0	EA
MS-46421	Disposable Standard Miller Blade, Size 1	EA
MS-46422	Disposable Standard Miller Blade, Size 2	EA
MS-46423	Disposable Standard Miller Blade, Size 3	EA
MS-46424	Disposable Standard Miller Blade, Size 4	EA
MS-466000	Reusable Standard Macintosh Blade, Size 00	EA
MS-46600	Reusable Standard Macintosh Blade, Size 0	EA
MS-46601	Reusable Standard Macintosh Blade, Size 1	EA
MS-46602	Reusable Standard Macintosh Blade, Size 2	EA
MS-46603	Reusable Standard Macintosh Blade, Size 3	EA
MS-46604	Reusable Standard Macintosh Blade, Size 4	EA
MS-466100	Reusable Standard Miller Blade, Size 00	EA
MS-46610	Reusable Standard Miller Blade, Size 0	EA
MS-46611	Reusable Standard Miller Blade, Size 1	EA
MS-46612	Reusable Standard Miller Blade, Size 2	EA
MS-46613	Reusable Standard Miller Blade, Size 3	EA
MS-46614	Reusable Standard Miller Blade, Size 4	EA
MS-46000	Handle, Standard, Adult	EA
MS-46001	Handle, Standard, Pediatric	EA


Features:

- Stainless Steel Blades
- Light Weight Design
- Smooth Surface Satin Finish Blades
- 2 x C Cell Batteries to Boost LED Light Output
- ISO 7376 Compatible:

Benefits:

- Max strength and anti-corrosive blades offer better performance and quality.
- Improves balance and maneuverability.
- Round smooth edges for greater patient comfort and satin finish reduces glare for airway visualization.
- LED bulbs offer a high quality light with a cool bulb for patient comfort.
- Complies with general requirements for laryngoscopes used for intubation, and specifies critical dimensions for the handle and lamp of hook-on type laryngoscopes.

Specification:



Brightest LED laryngoscope in the market



Strong Stainless Steel reinforcement in the base making it stronger than all plastic base laryngoscopes plastic jacket



Handles feature all metal heads.



Textured metal handles for maximum control and handling



2 "C-Cell" batteries included in handle to boost LED output

