



MedSource

I.V. Infusion

ProSite™ Male Luer PRN

MedSource PRN connectors are a reliable means of I.V. Administration. Needleless for patient and caregiver safety, MedSource is making safety affordable and safety compliance easy.

Our ProSite™ PRN is a Positive Displacement adapter. The PRN is made from Medical Grade PVC, DEHP / Not made with natural rubber latex and is Non-Pyrogenic.



Specification:

- 100/Case
- Priming Volume (mL) = 0.4 (max)
- Positive Displacement (mL) = .4 (max)
- Residual Volume (mL) = 0.35
- Flow Rate = 6.0 L per Hour.
- Negative Pressure Aspiration = No ingress of air up to 83 in -H20 of vacuum.
- Maximum System Pressure = No leakage up to 45 psi.
- Positive Pressure Flow at Disconnect = 0 mm to .030 mm
- Residual Volume Flushability through Access Port = 5 mL (max)
- Materials = Commonly used medical grade plastics compatible with sterilization processes; ISO 10993 certified.

Features:

- Needleless
- Non-Pyrogenic
- Latex/DEHP free
- Sterile, individually packaged, Swabable

Benefits:

- Minimizes potential for needle sticks.
- Will not cause inflammation.
- Safe for all patients.
- Prevents infection and cross contamination.

Rev. 1

www.gomedsource.com

Item #	Description	UOM
MS-83011	ProSite™ Male Luer PRN	EA, 100/CS