



**MedSource**  
Respiratory & Airway

# Aluminum O2 Cylinder

The MedSource Aluminum Cylinders are seamless and about 40% lighter weight than steel. They are made from high strength aluminum alloy 6061-T6.

There are two sizes to choose from and the standard finish is a brushed body with green crown and clear paint.

They are fully DOT approved and fully labeled.



## Specification:

	CAPACITY/Cu. ft/Liters	Length	Diameter	Empty Weight Lb/Kg	Thread Size
MS-MED-ET:	24.0/679	25.4/645	4.38/111.3	7.8/3.6	0.750 - 16UNF - 2B
MS-MED-DT:	15.0/425	16.7/424	4.38/111.3	5.4/2.4	0.750 - 16UNF - 2B

Rev. 1

[www.gomedsource.com](http://www.gomedsource.com)

Item #	Description	UOM
MS-MED-ET	Aluminum Oxygen Cylinder, E Tank	EA
MS-MED-DT	Aluminum Oxygen Cylinder, D Tank	EA