



MedSource

Transport &
Immobilization

Break-Apart Scoop Stretchers

Maximum Weight Capacity: 350lbs
Length: 67"- 80.5"
Width: 17"



MS-SCP123

Maximum Weight Capacity: 350lbs
Length: 63.5"-78"
Width: 18"



MS-SCP124

The MedSource Break-Apart Aluminum Scoop Stretcher is designed to be gently maneuvered under the patient without rolling or lifting. An open center permits the patient to be X-rayed while remaining immobilized on the stretcher.

Composed of lightweight aluminum with an adjustable length and 3 patient restraint straps.

The MedSource X-ray Translucent Break-Apart Scoop Stretcher is constructed of impact-resistant, lightweight, easy to clean polyethylene plastic. The head section is recessed to provide proper cervical alignment, allowing the patient's spine to remain in the position found.

Comes with 3 patient restraint straps.

Rev. 1

www.gomedsource.com

Item #	Description	UOM
MS-SCP123	Break-Apart Stretcher, Aluminum	EA
MS-SCP124	Break-Apart Stretcher, X-Ray Translucent	EA