



**MedSource**  
Wound Care

## EMT Shears



The large handle provides greater comfort while cutting. Safety tip prevents injury to the patient. Ideal for cutting through clothing, seatbelts, leather, soft metal, wire etc. Each shear measures 7 1/4 ". Available in 8 assorted colors.

### Specification:

- Hardened stainless steel
- Molded "comfort" handle
- Safety tip
- Autoclavable

Rev. 1

[www.gomedsource.com](http://www.gomedsource.com)

Item #	Description	UOM
MS-SH001B	EMT Shears, Black	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-SH002BL	EMT Shears, Blue	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-SH003P	EMT Shears, Hot Pink	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-SH004Y	EMT Shears, Yellow	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-SH005R	EMT Shears, Red	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-SH006G	EMT Shears, Green	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-SH007PR	EMT Shears, Purple	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-SH008O	EMT Shears, Orange	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS