



MedSource

Blankets and Linens

Cot Sheets



MS-001FL
Cot Sheet, Flat

- Durable, super strong sheet, 38 gsm
- Tear resistant
- Soft to the touch and comfortable for the patient



MS-001FS
**Cot Sheet, Fitted,
Sewn Ends**

- Made from a strong, durable emergency linen, 55 gsm
- Fluid impervious
- Specifically designed for emergency service stretchers
- Sewn corners for a secure fit



MS-42606
Cot Sheet, Fitted

- Made from a strong durable emergency linen with elastic corners for a secure fit
- Specifically designed for emergency service stretchers

Rev. 1

www.gomedsource.com

Item #	Description	UOM
MS-001FL	Cot Sheet, Flat, 87" x 40"	50, EA/Case
MS-001FS	Cot Sheet, Fitted, Sewn Ends, 72" x 6" x 30"	50, EA/Case
MS-42606	Cot Sheet, Fitted, Elastic Ends, 30" x 84"	50, EA/ Case



MedSource
Blankets and Linens

Heavy Duty Cot Sheets



**MS-42610
Cot Sheet, Heavy
Duty**

- Made from a strong, durable emergency linen, 45 gsm
- Provides fluid impervious barrier
- Weight tested to accommodate patient transfer



**MS-42604
Pro-fitted Cot Sheet,
Heavy Duty, Sewn
Ends**

- Made from a strong, durable emergency linen, 45 gsm
- Provides fluid impervious barrier
- Sewn corners for a secure fit



**MS-001HS
Cot Sheet, Heavy
Duty, Square Ends**

- Made from a spun-bonded polypropylene, 95gsm
- Fluid impervious
- Specifically designed for emergency service stretchers
- Elastic corners for a secure fit

Rev. 1

www.gomedsource.com

Item #	Description	UOM
MS-42610	Cot Sheet, Heavy Duty, Elastic Ends	25 EA, /Case
MS-42604	Heavy Cot Sheet Fitted Square Ends	50 EA, /Case
MS-001HS	Heavy Duty Cot Sheet, Fitted, Square Ends	50 EA, /Case



MedSource

Blankets and Linens

Speciality Cot Sheets



MS-44100 Poly Tissue Disposable Sheet

- Fluid resistant polyethylene outer surface provides strength and durability
- Tissue layer on the inside for comfort and softness.
- Lightweight yet durable



MS-42205 Heavy Duty Cot Sheet w/ Pouch

- Fluid impervious polyethylene, 60gsm
- Air laid paper under side for superior strength and support for the patient, 115gsm
- Elastic ends for secure fit
- Convenient pouch sewn into sheet



MS-42210 Ready Fit Fitted sheet

- Heavy Duty Cot Sheet specially designed for XPS Mattress.
- Fluid repellent for mattress protection
- Heavy duty construction for durability and strength
- Soft non-woven fabric for enhanced patient comfort

Rev. 1

www.gomedsource.com

Item #	Description	UOM
MS-44100	Poly-Tissue Flat Sheet, Economy	50 EA, /Case
MS-42205	Heavy Duty Cot Sheet with Pouch	36 EA, /Case
MS-42210	Ready Fit Cot Sheet	30EA, /Case