



**MedSource**  
Respiratory & Airway

# Robertazzi Nasopharyngeal Airway



Our Robertazzi Nasopharyngeal Airway is made from a soft, flexible neoprene material. The Latex-Free Robertazzi is the ideal way to maximize patient comfort while offering the best solution for latex sensitivity.

Rev. 1

[www.gomedsource.com](http://www.gomedsource.com)

Item	Description	UOM
MS-23950	Robertazzi Naso Airways, 20 Fr	10 EA/ BX
MS-23952	Robertazzi Naso Airways, 22 Fr	10 EA/ BX
MS-23954	Robertazzi Naso Airways, 24 Fr	10 EA/ BX
MS-23956	Robertazzi Naso Airways, 26 Fr	10 EA/ BX
MS-23958	Robertazzi Naso Airways, 28 Fr	10 EA/ BX
MS-23960	Robertazzi Naso Airways, 30 Fr	10 EA/ BX
MS-23962	Robertazzi Naso Airways, 32 Fr	10 EA/ BX
MS-23964	Robertazzi Naso Airways, 34 Fr	10 EA/ BX
MS-23966	Robertazzi Naso Airways, 36 Fr	10 EA/ BX

## Features:

- Softest Latex Free material
- Larger Trumpet
- Rounded, beveled tip
- Neoprene material construction

## Benefits:

- Helps prevent damage on the nasal passages
- Safe and secure placement
- Safe, gentle and comfortable insertion
- Strong and durable, resists collaps and kinking

## Specification:

- Available in 9 sizes: 20Fr- 36Fr
- This product is not made with natural rubber latex
- Sterile in unopened package
- Disposable, single use only