

# ChoiceMMed™ Fingertip Pulse Oximeter



The ChoiceMMed™ Fingertip Pulse Oximeter unit features a single directional OLED display with six display modes and adjustable brightness. The Pulse Oximeter also displays the SpO2, PR and Pulse Bar. The automatic power-off feature preserves battery life and is also compact and lightweight for on-the-go use.

Item #	Description	U of M
MS-74002-C	Fingertip Pulse Oximeter w/ SpO2, PR & Pulse Bar	EA

## Fingertip Pulse Oximeter

### Specifications

#### Features

- Dual OLED color displays, blood oxygen level, Pulse Rate and Pulse Bar
- Six display modes; adjustable brightness
- Low power consumption; automatic power-off
- 2 AAA batteries

- Operation Temperature: 41°F ~ 104°F
- Storage Temperature: -13°F ~ 158°F
- Operation Humidity: 15% ~ 93% no condensation
- Storage Humidity: ≤93% no condensation

#### Technical Specifications

SpO2 Measurement Range .....	70 ~ 100%
SpO2 Measurement Accuracy .....	70 ~ 100% ±2%; 0~ 69% no definition
PR Measurement Range .....	30 ~ 250 bpm
PR Measurement Accuracy .....	30 bpm ~ 99 bpm ±2%; 100 bpm ~ 250 bpm ±2%