



# TrueSafe® Safety IV Catheter




### RapidFlash™ Notched Needle

(Featured in sizes 20g - 24g)  
The notched needle promotes immediate confirmation of venipuncture helping to reduce hit and miss insertion attempts



### PTFE Catheter Material

High tensile strength catheter material



## SPECIFICATIONS

- PTFE catheter material
- Push button style needle retraction
- Notched needle (sizes 20g - 24g)
- Radiopaque
- Non-pyrogenic, not made with DEHP
- Not made with natural rubber latex
- Disposable
- Siliconized ground stainless steel needle
- Straight hub
- ISO compliant color coding
- Needle Stick Prevention and Safety Act (NSPA) of 2000 compliant



[www.medsourcelabs.com](http://www.medsourcelabs.com)

Item #	Description	UOM
MS-84414	TrueSafe® Safety IV Catheter, 14G x 1.75"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84416	TrueSafe® Safety IV Catheter, 16G x 1.77"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-844182	TrueSafe® Safety IV Catheter, 18G x 1.16"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84418	TrueSafe® Safety IV Catheter, 18G x 1.88"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-844202	TrueSafe® Safety IV Catheter, 20G x 1.16"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84420	TrueSafe® Safety IV Catheter, 20G x 1.88"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-844203	TrueSafe® Safety IV Catheter, 20G x 1.00"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84422	TrueSafe® Safety IV Catheter, 22G x 1.00"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84424	TrueSafe® Safety IV Catheter, 24G x 0.75"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS

# TrueSafe® Comfort Safety IV Catheter



**Triple Facet Back Cut Needle**  
Provides ease of insertion and reduced trauma to puncture site

**Push Off Tab**  
Assists with catheter advancement

**Push Button Safety Mechanism**  
Retracts needle fully and safely into catheter chamber

**Fully Encapsulating Needle Shield**  
Reduces exposure to blood borne pathogens via needle stick injuries

**Gently Tapered Catheter Tip**  
Provides smooth transition from needle to catheter during insertion

**Ridged Catheter Hub**  
Provides no slip grip during connection of IV Administration sets

**Contoured Housing**  
Ergonomically designed grips allow for one handed technique

**Large-Volume Chamber**  
Allows clear visual confirmation through flash-back and easy handling

**RapidFlash™ Notched Needle**  
(Featured in sizes 20g - 24g)  
The notched needle promotes immediate confirmation of venipuncture helping to reduce hit and miss insertion attempts

**VenaPUR™ Catheter Material**  
MedSource proprietary PUR material designed to remain kink free and softens within the body

## SPECIFICATIONS

- VenaPUR™ catheter material
- Push button style needle retraction
- Notched needle (sizes 20g - 24g)
- Radiopaque
- Non-pyrogenic, not made with DEHP
- Not made with natural rubber latex
- Disposable
- Siliconized ground stainless steel needle
- Straight hub
- ISO compliant color coding
- Needle Stick Prevention and Safety Act (NSPA) of 2000 compliant



www.medsourcelabs.com

Item #	Description	UOM
MS-84314	TrueSafe Comfort®, 14G x 1.75"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84316	TrueSafe Comfort®, 16G x 1.77"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-843182	TrueSafe Comfort®, 18G x 1.16"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84318	TrueSafe Comfort®, 18G x 1.88"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84320	TrueSafe Comfort®, 20G x 1.88"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-843202	TrueSafe Comfort®, 20G x 1.16"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-843203	TrueSafe Comfort®, 20G x 1.00"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84322	TrueSafe Comfort®, 22G x 1.00"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84324	TrueSafe Comfort®, 24G x 0.75"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-843242	TrueSafe Comfort®, 24G x 0.56"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84326	TrueSafe Comfort®, 26G x 0.56"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-843262	TrueSafe Comfort®, 26G x 0.75"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS