

Single Quart Sharps Container



Transparent Locking Lid

Lid locks to secure biohazard waste of needles, syringes and IV catheters

Translucent Red ABS Plastic With Fill level Identification

Puncture and leak resistant material is also transparent for fill level identification

Single Quart Size

Small size fits easily into gear bags, but large enough to use multiple times over

SPECIFICATION

Materials

ABS Plastic

Dimensions

Height: 6.25"

Depth: 4.25"

Width: 4.25"

Product Testing

Tested to the Standard of ISO ASTM F2132-01 for Puncture Resistance

ISO Standards

Meets ISO 23907:2012 Sharps Protection Standard

Shipping Carton Dimensions

30x32x54cm



www.medsourcelabs.com

Item #	Description	UOM
MS-64251	Single Quart Sharps Container	60 EA/CS