

# Latex-Free Tourniquet



**“Grippy” Surface**  
Better holds the patient’s  
arm when drawing blood

**Material Closely Resembles Latex**  
No worries about latex allergies

## SPECIFICATION

- 1” x 18”
- Materials Closely Resemble Latex
- Rolled or Unrolled

Item #	Description	UOM
MS-17301	Latex-Free Tourniquet, 1” x 18”, Disposable, Flat	1000 EA/CS
MS-17302	Latex-Free Tourniquet, 1” x 18”, Disposable, Rolled	1000 EA/CS