



SPECIFICATION

- MS-83091, MS-83092, MS-83094 Priming Volume: 2.0 mL
- MS-83097 Priming Volume: 3.2 mL
- MS-83216 Priming Volume: 2.3 mL
- Latex and DEHP Free
- Various Lengths and Configurations
- Manufactured with Medical Grade PVC

Item #	Description	UOM
MS-83090	Needleless Blue PRN, 6", Split Septum Y Site, Side Clamp	100/CS
MS-83091	Needleless Clear PRN, 6", Roberts Clamp, Rotary Male Luer Lock	100/CS
MS-83092	Needleless Blue PRN, 8", Roberts Clamp	100/CS
MS-83094	Needleless Blue PRN, 10", Roberts Clamp	100/CS
MS-83097	Flow Regulator, 19", Split Septum Y Site, Side Clamp	48/CS
MS-83216	Male Luer Lock, 21"	50/CS