

Cot Sheets


Special Use



MS-44100 Poly Tissue Disposable Sheet

SPECIFICATION

- Fluid resistant* polyethylene outer surface provides strength and durability
- Tissue on the inside for comfort and softness
- Lightweight yet durable



MS-42205 Heavy Duty, Cot Sheet with Pouch

SPECIFICATION

- Fluid resistant* polyethylene, 60 gsm
- Air laid underside for superior strength and support for the patient, 115 gsm
- Elastic ends for secure fit
- Convenient pouch sewn into sheet



MS-42210 Ready Fit, Fitted Sheet

SPECIFICATION

- Heavy Duty Cot Sheet made specially for XPS Mattress
- Fluid resistant* for mattress protection
- Heavy duty construction for durability and strength
- Soft non-woven fabric for enhanced patient comfort

**Fluid could pass through the material with pressure or time*



www.medsourcelabs.com

Item #	Description	UOM
MS-44100	Poly-Tissue Flat Sheet, Economy	50 EA/CS
MS-42205	Heavy Duty Cot Sheet with Pouch	36 EA/CS
MS-42210	Ready Fit Cot Sheet	30 EA/CS