




Our I.V Administration Sets are the first in the industry to introduce a color coded packaging system, to ensure you can always find the drip rate you are looking for.




SPECIFICATION:

- Individually packaged
- Sterile
- Not made with natural rubber latex
- DEHP Free
- Non-pyrogenic
- Medical Grade PVC

MS-83110 - 10 Drop
 MS-83115 - 15 Drop
 MS-83160 - 60 Drop


50/CASE



83", Non Vented Drip Chamber, Roller Clamp, Split Septum Y-Site, Needleless PRN, Rotary Male Luer Lock

MS-831106 - 10 Drop
 MS-831156 - 15 Drop
 MS-831606 - 60 Drop


50/CASE



83", Non Vented Drip Chamber, Roller Clamp, Split Septum Y-Site, Needleless PRN, Rotary Male Luer Lock
 With pre-attached 6" extension set

MS-83250 - 10 Drop


25/CASE



83", Vented Drip Chamber, Slide Clamp, Roller Clamp, Split Septum Y-Site, Needleless PRN, Stop Cock, Rotary Male Luer Lock

MS-83401 - 10 Drop
 MS-83406 - 15 Drop
 MS-83403 - 60 Drop


50/CASE



72", Non Vented Drip Chamber, Roller Clamp, Y-Site, Needleless PRN, Rotary Male Luer Lock

MS-83170T - 10 Drop
 MS-83171T - 60 Drop


50/CASE



88", Vented Drip Chamber, Roller Clamp, two (2) Roberts Clamps, Split Septum Y-Site, Needleless PRN, Rotary Luer Lock, Double Sided Tape

MS-83408 - 15 Drop

48/CASE



96", Non Vented Drip Chamber, Back Check Valve, Slide Clamp, Flow Regulator, two (2) Split Septum Y-Sites, Rotary Male Luer Lock


Item #	Description	UOM
MS-83110	83", 10 Drop, Split Septum Y-Site, Needleless PRN	50 EA/CS
MS-83115	83", 15 Drop, Split Septum Y-Site, Needleless PRN	50 EA/CS
MS-83160	83", 60 Drop, Split Septum Y-Site, Needleless PRN	50 EA/CS
MS-831106	83", 10 Drop, Split Septum Y-Site, Needleless PRN, 6" ext. set	50 EA/CS
MS-831156	83", 15 Drop, Split Septum Y-Site, Needleless PRN, 6" ext. set	50 EA/CS
MS-831606	83", 60 Drop, Split Septum Y-Site, Needleless PRN, 6" ext. set	50 EA/CS
MS-83250	83", 10 Drop, Split Septum Y-Site, Needleless PRN	25 EA/CS
MS-83401	72", 10 Drop, Y-Site, Needleless PRN	50 EA/CS
MS-83406	72", 15 Drop, Y-Site, Needleless PRN	50 EA/CS
MS-83403	72", 60 Drop, Y-Site, Needleless PRN	50 EA/CS
MS-83170T	88", 10 Drop, Split Septum Y-Site, Needleless PRN	50 EA/CS
MS-83171T	88", 60 Drop, Split Septum Y-Site, Needleless PRN	50 EA/CS
MS-83408	96", 15 Drop, 2 Split Septum Y-Sites	48 EA/CS

I.V. Administration Sets, 100+ inch

SPECIFICATION:

MS-83175

10 Drop
 Length: 100" (254cm)
 Not made with natural rubber latex.
 Non-pyrogenic,
 DEHP free,
 Single Use only
 Sterile in unopened, undamaged packaging




100", 10 Drop, Non Vented Drip Chamber,
 2 Slide Clamps, 1 Roller Clamp, 2 Split
 Septum Y-Sites and a Rotary Luer Lock

SPECIFICATION:

MS-83105

15 Drop
 Length: 105" (266cm)
 Not made with natural rubber latex.
 Non-pyrogenic,
 DEHP free,
 Single Use only
 Sterile in unopened, undamaged packaging




105", 15 Drop, Vented Drip Chamber,
 Back Check Valve, 1 Slide Clamp, 1 Roller
 Clamp, 2 Needleless PRN Y-Sites and
 a Rotary Luer Lock

SPECIFICATION:

MS-83122

15 Drop
 Length: 122" (309cm)
 Not made with natural rubber latex.
 Non-pyrogenic,
 DEHP free,
 Single Use only
 Sterile in unopened, undamaged packaging



122", 15 Drop, Vented Drip Chamber,
 Roller Clamp, Back Check Valve, Stop Cock,
 Needleless PRN Y-Sites and a Rotary
 Luer Lock

Item #	Description	U of M
MS-83175	IV Administration Solution Set, 100"	50 EA/Case
MS-83105	IV Administration Set 105", with back-check valve	50 EA/Case
MS-83122	IV Administration Set 122", with back-check valve and stopcock	50 EA/Case