

ChoiceMMed™ Fingertip Pulse Oximeter



SPECIFICATION

- SpO2 Measurement Range: 70 ~100%
 - SpO2 Measurement Accuracy: 70 ~ 100% $\pm 2\%$; 0 ~ 69% no definition
 - PR Measurement Range: 30 ~ 250 bpm
 - PR Measurement Accuracy: 30bpm ~ 99bpm $\pm 2\%$; 100bpm ~ 250bpm $\pm 2\%$
 - Dual OLED Color Displays, Blood Oxygen Level, Pulse Rate and Pulse Bar
 - Low Power Consumption; Automatic Power-Off
 - 2 AAA Batteries
-
- Operation Temperature: 41°F ~ 104°F
 - Storage Temperature: -13°F ~ 158°F
 - Operation Humidity: 15% ~ 93% no condensation
 - Storage Humidity: $\leq 93\%$ no condensation



www.medsourcelabs.com

Item #	Description	U of M
MS-74002-C	Fingertip Pulse Oximeter w/ SpO2, PR & Pulse Bar	EA